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- 1) Hardisty, R. M., and Macpherson, J.C: A one stage Factor VIII (antiheemophilic globulin) assay and its use on venous and capillary plasma. Thromb. Diath. haemorrh., 7.215, 1962
- 2) Douglas A. Triplett, MD & Cathy S. Harms, MT: Procedures for the Coagulation Laboratory. Educational Products Division American Society of Clinical Pathologists. P34-60.
- 3) 山田兼雄, 目黒 嵩: 血液凝固 - 止血と血栓・下. 宇宙堂 八木書店. P427~442

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## **FACTOR DEFICIENT PLASMAS (FII, FV, FVII, FVIII, FIX, FX, FXI, FXII)**

### **INTENDED USE:**

These plasmas are intended to be used as negative controls for coagulation factor assays or for preparing a standard curve for a coagulation factor assay. **FOR LABORATORY USE ONLY.** This plasma is from a human donor with the specific factor deficiency indicated on the vial. Different plasmas exhibit different sensitivities with different reagents. Each plasma should be evaluated for sensitivity with the system being used in the lab.

### **PREPARATION:**

Place vial in 37° C. water bath until plasma is thawed. Exact time is determined by volume of plasma in vial - approximately 2-5 minutes in circulating water bath. Mix gently and keep on ice until ready to use. Plasma should be used within 2 hours after being thawed.

### **WARNINGS:**

Individual donors used in this plasmas are tested negative for HBsAg, HIV and HCV by FDA approved test methods in the licensed blood centers where the plasmas were drawn. Because no test method can offer complete assurance that these or other infectious agents are absent, this product should be handled with universal precautions.

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